

JDGALAKIT00 CATARACT SET

5 PER BOX STERILE


★ BEST SELLER

set includes:

- . Ring handle scissors, straight, blunt
- . Adjustable speculum
- . Colibri forceps, 0.12 mm teeth
- . MICS Rhexis forceps, for incisions down to 1.8 mm, with controlled opening
- . Phaco chopper, aluminium handle



STERILE EO



JCP0003 TORIC IOL SET


5 PER BOX STERILE

set includes:

- . Ring + axial marke
- . Ring with marks each 5°



STERILE EO



JCPCAT1 CATARACT SET

5 PER BOX STERILE

set includes:


- . Ring scissors, straight, blunt tips
- . Adjustable speculum, 14 mm v-shaped blades
- . Corneal forceps, 1x2 teeth 0.20 mm
- . Capsulorrhexis forceps, angled, 11 mm, curved arms, sharp tips, ergonomic handle
- . Nagahara Phaco chopper blunt tip, 60° cutting edge, 10 mm angle to tip
- . Nucleus splitter-manipulator, delicate, 10 mm angle to tip



STERILE EO



5 Years Shelf Life



JCP0004.41 FLACS DISPOSABLE KIT

5 PER BOX STERILE


set includes:



STERILE EO



. Adjustable speculum



. Wire speculum




. Corneal forceps, 1x2 denti retti da 0.20 mm



. Incisions breaching spatula for femto laser procedure, 0.8 mm lancet shaped



. Tying forceps, angled, long, delicate, without platform



. Nucleus splitter-manipulator, delicate, 10 mm angle to tip



. IOL manipulator



. 90 mm scissors, straight, 24 mm long blades, pointed tips



. Viscoelastic cannula



. IOL lens loading forceps, steel



JDIOLEXP EXTRACTION IOL SET



5 PER BOX STERILE

set includes:

- . Scissors for removal of IOL, angled, 6 mm long, delicate blades, with serrated lower blade
- . IOL extraction forceps, angled, strongly serrated
- . Reverse hook



STERILE EO



JDFEMTO.1 FEMTO CATARACT BASIC SET

5 PER BOX STERILE

set includes:

- . Reverse speculum, solid blade. Medical polymer. Latex-free
- . Incisions breaching spatula for femto laser procedure, 0.8 mm lancet shaped (to open the incision performed by femto laser)
- . Nucleus splitter, disc shaped very thin for femto second laser (useful to separate the quadrants of the nucleus) disposable



STERILE EO




ADDITIONAL TOOLS

10 PER BOX STERILE



STERILE EO



JD2192.35A/10

Incisions breaching spatula for femto laser procedure, **0.8 mm** lancet shaped (to open the incision performed by femto laser)



JD2294.50/10

Anello graduato **5°** ogni step

Ring with marks each **5°**

ADDITIONAL TOOLS FOR THE BASIC SET JDFEMTO.1



STERILE EO



JD2180.24A/10

Nagahara Phaco chopper blunt tip, **60°** cutting edge, 10 mm angle to tip



JD3126.9/10

blunt tips

Capsule removal forceps with controlled opening, 11 mm **blunt tips**, work down to **1.7 mm** corneal incisions (to be used after femtosecond laser capsulotomy to grasp and remove the capsule)



STERILE EO




THIN-WALL



JD96624

23G infusion handpiece, "**thin-wall**", two **0.4 mm** sideports



JD96613

23G aspiration handpiece, curved, **sanblasted**, "**thin-wall**", **0.3 mm** port