

# JDLASIK01 LASIK SET


5 PER BOX STERILE

## set includes:

- . Hook to unstick the corneal flap. The delicate hook works at the flap edges inside the groove
- . Spatula to lift the corneal flap. The "C" shaped tapered spatula, once positioned under the flap, lifts it with a gentle push on the horizontal plane (from the hinge to the exit point)
- . Reverse speculum, solid blade. medical polymer. latex-free
- . LASIK irrigation cannula, flat, 3 ports, 25 gauge, 12 mm angle to tip
- . PVA spears with handle, 5 sponge



STERILE EO



# JCP0002 PRK PHOTO REFRACTIVE KERATOMY SET

5 PER BOX STERILE

## set includes:

- . Reverse speculum, solid blade. medical polymer, latex-free
- . Spatula for epithelium removal
- . Corneal marker, Ø 9.0 mm
- . PRK/LASIK irrigation cannula. Curved and flat, front opening



STERILE EO

