

DISPOSABLE MEIBOMIAN GLANDS PROBE

5 PER BOX STERILE




STERILE EO



5 Years
Shelf Life

DRY EYE

Probes for Meibomian glands
each size is provided in 5 units per box



SET

JD2123.2 probe length L=2mm - Ø=0.10mm

JD2123.4 probe length L=4mm - Ø=0.10mm

JD2123.6 probe length L=6mm - Ø=0.10mm

JD2124.2 probe length L=2mm - Ø=0.15mm

JD2124.4 probe length L=4mm - Ø=0.15mm

JD2124.6 probe length L=6mm - Ø=0.15mm